


PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

TRAVEL SCOOTERS





Contents

GO GO Mg-Lite	3
GO GO MG2	4
Apex Aluma-Lite 2.0	5
APEX Rapid	6
GO GO Ultra	7
GO GO Elite Traveller Sport	8
GO GO Elite Traveller 2	9
GO GO Endurance	10
GO GO Super Portable	11
Revo 2.0	12
GO GO Compact	13
Accessories	16
Comparison Table	18

Can hold up to
18 st 9 lbs

Can drive up to
4 mph

Weighs only
12 kg

Range up to
12.5 miles

Weighs
Only
12 kg

Digital display with
contactless lock

Battery level display

Front and rear
LED lighting



KEY FEATURES

- Lightweight quick release seat with adjustable armrests
- Easy disassembly
- Adjustable armrests
- Off board charging

Specifications

Weight capacity:	18 st 9 lbs (120 kg)
Overall length:	990 mm (38.9")
Overall width:	625 mm (24.6")
Front/Rear tyres:	Front - 195 mm (7.5") solid Rear - 195 mm (7.5") solid
Turning radius:	1450 mm (57")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	11 Ah Li (Flightsafe)
Battery range:	9.3 miles (15 km)
Total weight without batteries and seat:	12 kg
Battery weight:	1.6 kg each





Can hold up to
19 st 6 lbs



Can drive up to
4 mph



Weighs only
14.7 kg



Range up to
12.5 miles



Optional storage
available



Lightweight lithium
batteries



Foldable compact
design



Can hold up to
18 st 9 lbs



Can drive up to
4 mph



Weighs only
31.8 kg



Range up to
11.5 miles



New delta tiller



Airline safe batteries



8" alloy wheels

GOGO M2

**Motability
Scheme**

- Optional arm rests available
- 2 step easy fold system
- Optional battery sizes available
- Available in two different colours
- Compact lightweight design

Specifications

Weight capacity:	19 st 6 lbs (124 kg)
Overall length:	1016 mm (40")
Overall width:	429 mm (16.9")
Front/Rear tyres:	Front - 203 mm (8") solid Rear - 203 mm (8") solid
Turning radius:	1200 mm (47")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	6 Ah Li: (optional) 10 Ah Li: (standard) 20 Ah Li: (optional)
Battery range:	6 Ah Li: 4 miles (6.4 km) 10 Ah Li: 6.2 miles (10 km) 20 Ah Li: 12.5 miles (20 km)
Total weight without batteries:	14.7 kg
Battery weight:	6 Ah Li: 1 kg 10 Ah Li: 1.7 kg 20 Ah Li: 3 kg



APEX
ALUMA-LITE 2.0

**Motability
Scheme**

KEY FEATURES

- New extended footplate for extra leg room

Specifications

Weight capacity:	18 st 9 lbs (120 kg)
Overall length:	1022 mm (40.2")
Overall width:	490 mm (19.29")
Front/Rear tyres:	Front - 203 mm (8") solid Rear - 203 mm (8") solid
Turning radius:	1385 mm (54.5")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	12 Ah Li (flightsafe)
Battery range:	11.5 miles (18.5 km)
Total weight without batteries:	31.8 kg
Battery weight:	2.5 kg





Can hold up to
21 st 6 lbs



Can drive up to
4 mph




Weighs only
41 kg



Range up to
10.5 miles



Colour range



LED headlight



Ergonomic delta tiller



Can hold up to
18 st 9 lbs



Can drive up to
4 mph



Weighs only
36.3 kg



Range up to
11.5 miles



New longer footplate
for additional leg room



New delta tiller



8" alloy wheels



KEY FEATURES

- Disassembles into 5 pieces
- Front and rear suspension
- Optional 18 Ah batteries
- 8" alloy style wheels
- Compact fold flat seat

Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	1060 mm (41.75")
Overall width:	521 mm (20.5")
Front/Rear tyres:	Front - 203 mm (8") solid Rear - 203 mm (8") solid
Turning radius:	1219 mm (48")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 12 Ah: (standard) 2 x 18 Ah: (optional)
Battery range:	12 Ah: 5.5 miles (8.75 km) 18 Ah: 10.5 miles (17 km)
Total weight without batteries:	41 kg
Battery Weight:	12 Ah: 9 kg 18 Ah: 14 kg



(Standard colours)



KEY FEATURES

- 2 step easy fold system
- Optional battery sizes available
- Available in two different colours
- Budget friendly

Specifications

Weight capacity:	18 st 9 lbs (120 kg)
Overall length:	1022 mm (40.23")
Overall width:	490 mm (19.29")
Front/Rear tyres:	Front - 203 mm (8") solid Rear - 203 mm (8") solid
Turning radius:	1385 mm (54.5")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 12 Ah sealed lead-acid
Battery range:	11.5 miles (18.5 km)
Total weight without batteries:	36.3 kg
Battery weight:	4 kg





Can hold up to
23 st 3 lbs



Can drive up to
4 mph



Weighs only
44 kg



Range up to
14 miles



Feather-touch
disassembly



Front and rear
LED lighting



Extra 2" footplate
for extra leg room



Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
37.7 kg



Range up to
14 miles

Our
smallest
turning
radius
ever!



37" turning radius



Easy assembly



Storage comes
as standard



KEY FEATURES

- Patented EZ Technology offering a 37" turning radius
- Wraparound delta tiller
- 9" chunky solid tyres
- Off-board battery charging
- Compact fold flat seat



Specifications

Weight capacity:	23 st 3 lbs (147 kg)
Overall length:	1092 mm (43")
Overall width:	540 mm (21.25")
Front/Rear tyres:	Front - 229 mm (9") solid Rear - 229 mm (9") solid
Turning radius:	940 mm (37")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 20 Ah
Battery range:	20 Ah: 14 miles (23 km)
Total weight without batteries:	44 kg
Battery weight:	12.5 kg




KEY FEATURES

- Patented EZ Technology offering 37" turning radius
- New stylish mesh seat for extra comfort
- Arm rest storage bag comes as standard



Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	1000 - 1036 mm (39.37" - 40.79")
Overall width:	494 mm (19.45")
Front/Rear tyres:	Front - 178 mm (7") solid Rear - 203 mm (8") solid
Turning radius:	940 mm (37")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 12 Ah 2 x 18 Ah (optional)
Battery range:	12 Ah: Up to 11 miles 18 Ah: Up to 14 miles
Total weight without batteries:	37.7 kg (83 lbs)
Battery Weight:	12 Ah: 8.75 kg (standard) 18 Ah: 13.4 kg (optional)





Can hold up to
23 st 3 lbs



Can drive up to
4 mph



Weighs only
41 kg




Range up to
13.2 miles



Standard 8 Ah
airline-safe
Lithium battery



Dual charging ports
for both on and off
board charging



Delta tiller with
ergonomic wraparound
handles



Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
18.9 kg



Range up to
10 miles




360 degrees magnetic
charging point



Rear footpeg for upright
storage option



Storage bag
as standard



KEY FEATURES

- Compact, foldable seat with armrest
- Delta tiller controls
- High-Intensity LED headlight
- Interchangeable colour panels
- CTS suspension
- Feather-touch disassembly



Long lasting battery life
with Lithium batteries

Specifications

Weight capacity:	23 st 3 lbs (147 kg)
Overall length:	1070 mm (42")
Overall width:	690 mm (27")
Front/Rear tyres:	Front - 203 mm (8") solid Rear - 203 mm (8") solid
Turning radius:	1225 mm (48.25")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	8 Ah: (standard) 16 Ah (optional)
Battery range:	8 Ah: 6.6 miles (8 km) 16 Ah: 13.2 miles (21.3 km)
Total weight without batteries:	41 kg
Battery Weight:	8 Ah: 3.5 kg 16 Ah: 5.5 kg



KEY FEATURES

- Underseat storage
- Optional arm rests available
- Front and rear lighting
- 360 degree magnetic charging port
- Adjustable front tiller

Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	945 mm (37.2")
Overall width:	547 mm (21.54")
Front/Rear tyres:	Front - 178 mm (7") solid Rear - 203 mm (8") solid
Turning radius:	1311 mm (51.62")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	11.7 Ah Li (standard) flightsafe 16 Ah Li (optional)
Battery range:	Up to 10 miles (16.3 km)
Total weight without batteries:	18.9 kg
Battery weight:	11.7 Ah: 3.6 kg 16 Ah: 4.4 kg





Can hold up to
28 st 8 lbs




Can drive up to
5 mph




Weighs only
60 kg



Range up to
17 miles



Comfort trac front and
rear suspension



Underseat storage



LED headlight



Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
19 kg



Range up to
19 miles



Contactless locking
LED display



Anti theft alarm



Stay connected with the
Smart Mobility App

REVO 2.0



KEY FEATURES

- Convenient tiller charging point
- Easily disassembled into five pieces
- Storage compartment
- Monoshock rear suspension
- Available with either solid or pneumatic tyres
- Fold-flat seat with adjustable armrests

Specifications

Weight capacity:	28 st 8 lbs (181 kg)
Overall length:	1194 mm (47")
Overall width:	525 mm (20.75")
Front/Rear tyres:	Front - 229 mm (9") solid or pneumatic Rear - 229 mm (9") solid or pneumatic
Turning radius:	1323 mm (52")
Max speed:	5 mph (8 km/h)
Battery requirement:	30-35 Ah
Battery range:	Solid: Up to 17 miles (28.6 km) pneumatic: Up to 16.2 miles (26.1 km)
Total weight without batteries:	60 kg
Battery weight:	11.25 kg



Solid wheels available
in blue only

GO GO COMPACT




SCAN ME

KEY FEATURES

- Underseat storage
- Foldaway armrest
- Full Lighting package
- USB device charging port
- Cruise control



Innovative folding mechanism



Enhance your ride with our innovative app-designed for your mobility scooter.

Check your battery life, personalise your settings, and access useful features, all from your smartphone.

Mobility just got smarter!


Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	998 mm (39.29")
Overall width:	565 mm (22.24")
Front/Rear tyres:	Front - 203 mm (8") solid Rear - 229 mm (9") solid
Turning Radius:	1380 mm (54.3")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	10 Ah Li (Standard) Flightsafe 20 Ah Li: (Optional)
Battery range:	10 Ah Li: Up to 9.3 miles (15 km) 20 Ah Li: Up to 19 miles (30 km)
Total weight without batteries:	19 kg
Battery weight:	10 Ah: 3 kg 20 Ah: 6 kg





Stay connected with the
Smart Mobility App



UNLOCK THE FEATURES

1. Lock and unlock button to turn on and off the P mode. (Parking Mode)
2. Light button to turn on and off the front and tail lights.
3. Cruise Mode button to turn on and off the cruise control.

And much more...



Accessories



Rear view mirror

Optional rear view mirror for enhanced visibility and safety.



Foldable travel bags

Light and durable travel bags for the Quest, iGo, iGo Plus and iGo Fold.



Baskets & bags

Our handy and versatile range of additional large baskets and storage bags make transporting your essentials hassle free.



Crutch & Walker Holders

Quick and easy access for your cane, crutch or walking stick with our selection of holders.

Accessories



Rear oxygen tank holder

Featuring 2 thumb screws for safe transportation.



Replacement keys

We offer a replacement key service.



Cup holders

Satisfy your thirst on the move with our handy drinks holders, compatible with most of our range of products.



XLR USB charger & phone mounts

Charge your phone on the go!
This can be used while your powerchair or scooter is in motion.

Comparison Table

*All specifications may vary slightly. For more infomation on all our products, please check the product specification sheet found online.

	Weight Capacity	Overall Length	Overall Width	Front Tyres	Back Tyres	Turning Radius	Max Speed
	18 st 9 lbs (120 kg)	990 mm (38.9")	625 mm (24.6")	195 mm (7.5") solid	195 mm (7.5") solid	1450 mm (57")	4 mph (6.4 km/h)
	19 st 6 lbs (124 kg)	1016 mm (40")	429 mm (16.9")	203 mm (8") solid	203 mm (8") solid	1200 mm (47")	4 mph (6.4 km/h)
	18 st 9 lbs (120 kg)	1022 mm (40.2")	490 mm (19.29")	203 mm (8") solid	203 mm (8") solid	1385 mm (54.5")	4 mph (6.4 km/h)
	21 st 6 lbs (136 kg)	1060 mm (41.75")	521 mm (20.5")	203 mm (8") solid	203 mm (8") solid	1219 mm (48")	4 mph (6.4 km/h)
	18 st 9 lbs (120 kg)	1022 mm (40.2")	490 mm (19.29")	203 mm (8") solid	203 mm (8") solid	1385 mm (54.5")	4 mph (6.4 km/h)
	23 st 3 lbs (147 kg)	1092 mm (43")	540 mm (21.25")	229 mm (9") solid	229 mm (9") solid	940 mm (37")	4 mph (6.4 km/h)
	21 stone 6 lbs (136 kg)	1000 - 1036 mm (39.37" - 40.79")	494 mm (19.45")	178 mm (7") solid	203 mm (8") solid	940 mm (37")	4 mph (6.4 km/h)
	23 st 3 lbs (147 kg)	1070 mm (42")	690 mm (27")	203 mm solid (8")	203 mm (8")	1225 mm (48.25")	4 mph (6.4 km/h)
	21 st 6 lbs (136 kg)	945 mm (37.2")	547 mm (21.54")	178 mm solid (7")	203 mm solid (8")	1311 mm (51.62")	4 mph (6.4 km/h)
	28 st 8 lbs (181 kg)	1194 mm (47")	525 mm (20.75")	229 mm (9") solid or pneumatic	229 mm (9") solid or pneumatic	1323 mm (52")	5 mph (8 km/h)
	21 st 6 lbs (136 kg)	998 mm (39.29")	565 mm (22.24")	203 mm (8") solid	229 mm (9") solid	1380 mm (54.3")	4 mph (6.4 km/h)


Battery Requirements	Battery Range	Total weight Without Batteries	Battery Weight
12 Ah sealed lead acid	11.5 miles (18.5 km)	36.3 kg	4 kg
6 Ah Li: (Optional) 10 Ah Li: (Standard) 20 Ah Li: (Optional)	6 Ah Li: 4 miles 10 Ah Li: 6.2 miles (10 km) 20 Ah Li: 12.42 miles (20 km)	14.7 kg	6 Ah Li: 1 kg 10 Ah Li: 1.7 kg 20 Ah Li: 3 kg
12 Ah Li (flightsafe)	11.5 miles (18.5 km)	31.8 kg	2.5 kg
2 x 12 Ah (standard) 2 x 18 Ah: (optional)	2 x 12 Ah 5.5 miles (8.75 km) 18 Ah: 10.5 miles (17 km)	41 kg	12 Ah: 9 kg 18 Ah: 14 kg
2 x 12 Ah sealed lead-acid	11.5 miles (18.5 km)	36.3 kg	4 kg
2 x 20 Ah	20 Ah: 14 miles (23 km)	44 kg	12.5 kg
2 x 12 Ah 2 x 18 Ah (optional)	12 Ah: Up to 11 miles 18 Ah: Up to 14 miles	37.7 kg	12 Ah: 8.75 kg (standard) 18 Ah: 13.4 kg (optional)
12 Ah Li (Flightsafe)	11.5 miles (18.5 km)	31.8 kg	2.5 kg
11.7 Ah Li (standard) flightsafe 16 Ah Li (optional)	Up to 10 miles (16.3 km)	18.9 kg	11.7 Ah: 3.6 kg 16 Ah: 4.4 kg
30-45 Ah	Solid: Up to 17 miles (28.6 km) pneumatic: Up to 16.2 miles (26.1 km)	60 kg	11.25 kg
10 Ah Li: (standard) flightsafe 20 Ah Li: (optional)	10 Ah Li: Up to 9.3 miles (15 km) 20 Ah Li: Up to 19 miles (30 km)	19 kg	10 Ah: 3 kg 20 Ah: 6 kg

PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

01869 324600 sales@pride-mobility.co.uk pride-mobility.co.uk



Large Scooter Range




Powered Wheelchairs Range



Sitting Pretty
Rise and Recline chairs



Sitting Pretty
Adjustable Bed Range



Viva Sleep
Bed Range



REV 12 - 01/04/2026



Pride Mobility Products LTD 32 Wedgwood Road
Bicester, Oxon OX26 4UL ENGLAND