

PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

POWERED WHEELCHAIRS





Contents

Jazzy Carbon Ultra-Lite	3
Jazzy Carbon Swift	4
Jazzy Carbon Injection	5
Jazzy Carbon Max	6
Jazzy Carbon Max HD	7
Jazzy Select	8
Jazzy Select 6	9
Jazzy 600 ES	10
Jazzy Air 2	11
Fusion e	12
Go Chair	13
Comparison Table	14



Can hold up to
21 st 5 lbs



Can drive up to
4 mph



Weights only
14.5 kg



Range up to
18 miles

Weights
only
14.5 kg



One touch
fold function



Easy to use controls



Fold away footplate



- Mesh Seat for extra comfort
- Weighs only 14.5kg

Specifications

Weight capacity:	21 st 5 lbs (138 kg)
Overall length:	1000 mm (39.37")
Overall width:	620 mm (24.4")
Front/Rear tyres:	Front: 203 mm (8") solid Rear: 304 mm (12") solid
Turning radius:	900 mm (35.4")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	10 Ah Li (standard) 20 Ah Li (optional)
Battery range:	10 Ah: Up to 9 miles (14.8 km) 20 Ah: Up to 18 miles (28.96 km)
Total weight without batteries:	14.5 kg
Battery weight:	10 Ah: 1.4 kg (standard) 20 Ah: 3 kg (optional)



Can hold up to
18 st 10 lbs



Can drive up to
4 mph



Weighs only
10.26 kg



Range up to
9.3 miles

Weighs
only
10.26 kg



Narrow frame
offering effective
maneuverability



Flight safe batteries

Easy-fold mechanism



Can hold up to
21 st 4 lbs



Can drive up to
4 mph



Weighs only
19 kg



Range up to
9.3 miles



Flight safe batteries



Easy-to-use controls

Jazzy
Carbon

SWIFT



KEY FEATURES

- Easy-to-use controls
- Compact design

Specifications

Weight capacity:	18 st 10 lbs (120 kg)
Overall length:	800 mm (31.4")
Overall width:	565 mm (22.24")
Front/Rear tyres:	Front: 125 mm (5") solid Rear: 165 mm (6.5") solid
Turning radius:	900 mm (35.4")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	10 Ah Lithium (flightsafe)
Battery range:	9.3 miles (15 km)
Total weight without batteries:	10.26 kg
Battery weight:	1.6 kg

Jazzy
Carbon

INJECTION



KEY FEATURES

- Easy-fold mechanism
- New and improved cushion
- Honeycomb injection for improved strength

Specifications

Weight capacity:	21 st 4 lbs (136 kg)
Overall length:	950 mm (37.4")
Overall width:	610 mm (24")
Front/Rear tyres:	Front: 178 mm (7") solid Rear: 216 mm (8.5") solid
Turning radius:	900 mm (35.4")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	12 Ah Lithium (flightsafe)
Battery range:	9.3 miles (15 km)
Total weight without batteries:	19 kg
Battery weight:	2.1 kg



Can hold up to
28 st 7 lbs



Can drive up to
4 mph



Weighs only
12.6 kg



Range up to
9.3 miles



Heavy duty weight
capacity of up to 28.7 st



Lightweight strong
carbon fibre frame



Folding footplate



Can hold up to
28 st 6 lbs



Can drive up to
4 mph



Weighs only
13.8 kg



Range up to
12.4 miles

Heavy
duty
wider
seat



Heavy duty
wider seat



Lightweight,
pavement friendly



Folding footplate

Jazzy
CARBON MAX

**Motability
Scheme**

KEY FEATURES

- Easy folding mechanism
- Heavy duty weight capacity

Specifications

Weight capacity:	28 st 7 lbs (182 kg)
Overall length:	930 mm (36.6")
Overall width:	595 mm (23.4")
Front/Rear tyres:	Front: 175 mm (7") solid Rear: 210 mm (8) solid
Max speed:	4 mph (6.4 km/h)
Battery requirement:	10 Ah Li (Standard) Flightsafe 20 Ah Li: (Optional)
Battery range:	9.3 miles
Total weight without batteries:	12.6 kg
Battery weight:	2.1 kg



Jazzy
CARBON MAX HD

**Motability
Scheme**

KEY FEATURES

- Easy folding mechanism
- Heavy duty weight capacity

Specifications

Weight capacity:	28 st 6 lbs (182 kg)
Overall length:	965 mm (37.9")
Overall width:	665 mm (26.1")
Seat width:	510 mm x 470 mm
Front/Rear tyres:	Front: 175 mm (7") solid Rear: 270 mm (10") solid
Max speed:	4 mph (6.4 km/h)
Battery requirement:	24V 20Ah (DMEGC Cell)
Battery range:	12.4 miles
Total weight without batteries:	13.8 kg
Battery weight:	3.4 kg





Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
42.5 kg



Range up to
19.4 km



Active-Trac®
suspension



Seat depth
adjustable



Flip out footplate



Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
44.5 kg



Range up to
13 miles



Foldaway footplate



ATX Suspension
(Active-Trac With
Extra Stability)



Easy to use controller



KEY FEATURES

- Intuitive easy to operate controller
- Chunky 10" drive tyres

Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	1118 mm (44")
Overall width:	603 mm (23.75")
Front/Rear tyres:	Front: 146 mm (5.75") solid Rear: 254 mm (10") solid
Turning radius:	571 mm (22.5")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 12 volt
Battery range:	Up to 19.4 miles (31.2 km)
Total weight without batteries:	42.5 kg
Battery weight:	11.5 kg



KEY FEATURES

- Chunky 10" drive tyres
- Extra durable shrouds

Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	1067 mm (42")
Overall width:	597 mm (23.8")
Front/Rear tyres:	Front: 146 mm (5.75") solid Rear: 254 mm (10") solid
Turning radius:	546 mm (21.5")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 12 Volt deep cycle
Battery range:	Up to 13 miles (21 km)
Total weight without batteries:	44.25 kg
Battery Weight:	11.5 kg





Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
56 kg



Range up to
18.64 miles



Fold flat seat



Extra storage



Angle, height and
depth adjustable
armrest



Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
91.5 kg



Range up to
16.6 miles



Drives safely at
4 mph while
elevated at 12"



Fold flat seat



Angle and height
adjustable headrest

Jazzy
600 ES



KEY FEATURES

- Exceptional handling over tough terrain
- OMNI - casters on front and rear wheels to prevent hang-ups
- Active-Trac® suspension

Seat sizes available:

1. 16"W - 16" - 18"
2. 18"W - 18" - 20"
3. 20"W - 18" - 20"



Crash
tested

Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	1022 mm (40.25")
Overall width:	622 mm (24.5")
Front/Rear tyres:	Front: 152 mm (6") solid Rear: 152 mm (6") solid
Turning radius:	521 mm (20.5")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 50 Ah
Battery range:	Up to 18.64 miles (30 km)
Total weight without batteries:	56 kg
Battery weight:	17 kg



Jazzy
air₂



KEY FEATURES


- Storage compartment
- Height adjustable footrest
- Bright LED marker lights
- Active-Trac® suspension with mid-wheel 6® technology
- High back reclining sport seat

Seat sizes available:

1. 16"W - 16" - 18"
2. 18"W - 18" - 20"
3. 20"W - 18" - 20"

Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	1067 mm (42")
Overall width:	590 mm (23.25")
Front/Rear tyres:	Front: 152 mm (6") solid Rear: 254 mm (10") solid
Turning radius:	616 mm (24.25")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 34 Ah
Battery range:	Up to 16.6 miles (26.7 km)
Total weight without batteries:	91.5 kg
Battery Weight:	14.5 kg





Can hold up to
23 st 6 lbs



Can drive up to
4 mph



Weighs only
86 kg



Range up to
15.7 miles



Rear suspension



Optional lights
and indicators



VR2 control system



Can hold up to
21 st 6 lbs



Can drive up to
4 mph



Weighs only
44.25 kg



Range up to
13.2 miles



Disassembles into
five pieces



Large comfortable
footplate



Adjustable height
and width armrest

FUSION^e



KEY FEATURES

- Solid or Pneumatic tyres
- Independent rear wheel drive system
- Seat recline

Specifications

Weight capacity:	23 st 6 lbs (150 kg)
Overall length:	1050 mm (41.33")
Overall width:	630 mm (24.8")
Front/Rear tyres:	Front: 254 mm (10") solid Rear: 355 mm (14") solid
Turning radius:	641.5 mm (25.25")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 50 Ah
Battery range:	Up to 15.7 miles (25.2 km)
Total weight without batteries:	86 kg
Battery weight:	17 kg

GOCHAIR




KEY FEATURES

- Height and width adjustable armrests
- Dual swing-away under seat storage
- Large footplate


Specifications

Weight capacity:	21 st 6 lbs (136 kg)
Overall length:	842 mm (33")
Overall width:	575 mm (22.5")
Front/Rear tyres:	Front: 152 mm (6") solid Rear: 229 mm (9") solid
Turning radius:	648 mm (25.25")
Max speed:	4 mph (6.4 km/h)
Battery requirement:	2 x 18 Ah
Battery range:	18 Ah: 13.2 miles (21.3 km)
Total weight without batteries:	44.25 kg
Battery Weight:	14 kg



Comparison Table

*All specifications may vary slightly. For more infomation on all our products, please check the product specification sheet found online.

	Weight Capacity	Overall Length	Overall Width	Front Tyres	Back Tyres	Turning Radius	Max Speed
	21 st 5 lbs (138 kg)	1000 mm (39.37")	620 mm (24.4")	203 mm (8") solid	304 mm (12") solid	900 mm (35.4")	4 mph (6.4 km/h)
	18 st 10 lbs (120 kg)	800 mm (31.4")	565 mm (22.24")	125 mm (5") solid	165 mm (6.5") solid	900 mm (35.4")	4 mph (6.4 km/h)
	21 st 4 lbs (136 kg)	950 mm (37.4")	610 mm (24")	178 mm (7") solid	216 mm (8.5") solid	900 mm (35.4")	4 mph (6.4 km/h)
	28 st 7 lbs (182 kg)	930 mm (36.6")	595 mm (23.4")	175 mm (7") solid	210 mm (8) solid	850 mm (25.5")	4 mph (6.4 km/h)
	28 st 6 lbs (182 kg)	965 mm (37.9")	665 mm (26.1")	175 mm (7") solid	270 mm (10") solid	925 mm (36.4")	4 mph (6.4 km/h)
	21 st 6 lbs (136 kg)	1118 mm (44")	603 mm (23.75")	146 mm (5.75") solid	254 mm (10") solid	571 mm (22.5")	4 mph (6.4 km/h)
	21 st 6 lbs (136 kg)	1067 mm (42")	597 mm (23.8")	146 mm (5.75") solid	254 mm (10") solid	546 mm (21.5")	4 mph (6.4 km/h)
	21 st 6 lbs (136 kg)	1022 mm (40.25")	622 mm (24.5")	152 mm (6") solid	152 mm (6") solid	521 mm (20.5")	4 mph (6.4 km/h)
	21 st 6 lbs (136 kg)	1067 mm (42")	590 mm (23.25")	152 mm (6") solid	254 mm (10") solid	616 mm (24.25")	4 mph (6.4 km/h)
	23 st 6 lbs (150 kg)	1050 mm (41.33")	630 mm (24.8")	254 mm (10") solid	355 mm (14") solid	641.5 mm (25.25")	4 mph (6.4 km/h)
	21 st 6 lbs (136 kg)	842 mm (33")	575 mm (22.5")	152 mm (6") solid	229 mm (9") solid	648 mm (25.5")	4 mph (6.4 km/h)

Battery Requirements	Battery Range	Total weight Without Batteries	Battery Weight
10 Ah Li (standard) 20 Ah Li (optional)	Up to 18 miles (28.96 km)	14.5 kg	10 Ah: 1.4 kg (standard) 20 Ah: 3 kg (optional)
10 Ah Lithium (flightsafe)	9.3 miles (15 km)	10.26 kg	1.6 kg
12 Ah Lithium (flightsafe)	9.3 miles (15 km)	19 kg	2.1 kg
12 Ah Lithium (flightsafe)	9.3 miles (15 km)	12.6 kg	2.1 kg
24V 20Ah (DMEGC Cell)	12.4 miles (20 km)	13.8 kg	3.4 kg
2 x 12 volt	Up to 19.4 miles (31.2 km)	42.5 kg (base only)	11.5 kg
2 x 12 Volt	Up to 13 miles (21 km)	44.25 kg (base only)	11.5 kg
2 x 50 Ah	Up to 18.64 miles (30 km)	56 kg (base only)	17 kg
2 x 34 Ah	Up to 16.6 miles (26.7 km)	91.5 kg	14.5 kg
2 x 50 Ah	Up to 15.7 miles (25.2 km)	86 kg	17 kg
2 x 18 Ah	Up to 13.2 miles (21.3 km)	44.25 kg	14 kg

PRIDE[®]


MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

01869 324600 sales@pride-mobility.co.uk pride-mobility.co.uk



Large Scooter Range




Travel Scooter Range



Sitting Pretty
Rise and Recline chairs



Sitting Pretty
Adjustable Bed Range



Viva Sleep
Bed Range



REV 11 - 28/01/2026



Pride Mobility Products LTD 32 Wedgwood Road
Bicester, Oxon OX26 4UL ENGLAND